Easy Plate TM

Save time, labor, and cost with Kikkoman Easy Plate. Easy Plate is a convenient, prepared media plate for faster and easier microbiological testing.





Product Incubation Shelf life

Easy Plate SA lococcus aureus 35 or 37±1°C, 24±1 hours 12 months after manufacturing date Certification AOAC PTM, MicroVal, NordVal



Shelf life

Yeast and Mold rapid 25°C±1, 48±2 hours 18 months after manufacturing date



24±1 hours

12 months after

manufacturing date

Incubation

Shelf life

Colony Counting System for Easy Plate

Software Included

Automatically counts the colonies in imaged Easy Plate. Simpler and time-saving way to count colonies than manual counting advances efficiency and eliminates the effects of variation and skill level among operators.

Boost lab productivity easily with automatic colony counter



- Up to 7 sheets can be inserted at once.
- Highly accurate analysis with Al technologies.
- Automatically determines which of the 5 Easy Plate products is used and measures colonies.
- It also identifies spread colonies.
- Parameters to judge colonies are available with no initial settings.
- Able to modify the result and add the comments.

Supported Languages:

English, German, Spanish, French, Korean, Thai, Vietnamese, Simplified Chinese, Traditional Chinese, Japanese



*It requires the scanner (ADS-4300N by Brother Industries, LTD.) for imaging the Easy Plate. Available for purchase.



Kikkoman Biochemifa Company

2-1-1 Nishi-shinbashi, Minato-ku, Tokyo 105-0003 JAPAN biochemifa@mail.kikkoman.co.jp https://biochemifa.kikkoman.com/e/

© 2023 Kikkoman Corp. (6012W240101)





FMCG Industry Solutions Pty Ltd ABN 22 135 446 007

- Factory 11A, 1 3 Endeavour Rd, Caringbah NSW 2229
- www.fmcgis.com.au
- sales@fmcgis.com.au
- 얗 1300 628 104 or 02 9540 2288