

# PRO-Clean™

## Rapid Protein Residue Test

PRO-Clean is a quick and easy way to accurately monitor the cleanliness of surfaces to help ensure product quality. PRO-Clean detects protein residues left on a surface after cleaning. Simply swab a surface, release the reagent and if protein residue is present the reagent will turn purple. The color change provides a semiquantitative measure of the surface cleanliness. The more contamination present, the quicker the color change to purple and the darker the color. PRO-Clean quickly validates the hygiene of a surface, allowing immediate corrective action to be taken if necessary.



## DETECTS PROTEIN RESIDUES

### Key Features

- ✓ Results in under 10 minutes
- ✓ All-in-one device — no instrumentation
- ✓ 18-month shelf life at refrigerated temperatures (2 - 25°C)
- ✓ Easy to interpret
- ✓ Affordable


**NO INSTRUMENTATION**

### Ideal for

- ✓ Meat processors
- ✓ Food processors
- ✓ Restaurants and food service
- ✓ Grocery stores
- ✓ Industrial applications

**RESULTS IN**

**LESS THAN 10 MIN**



FMCG Industry Solution Ptv Ltd

ABN 22 135 446 007

A: Unit 11A, 1 - 3 Endeavour Rd, Caringbah NSW 2229

W: [www.fmcgis.com.au](http://www.fmcgis.com.au) E: [sales@fmcgis.com.au](mailto:sales@fmcgis.com.au)

Phone: 1300 628 104 / 02 9540 2288

ORDER # PRO-100

(100 tests per box)